



Code: SG1 – Seatgard

## quick start instructions

**Ramblegard**<sup>®</sup>  
stay a step ahead



1. Remove the monitor from the protective holder.
2. Open the Battery Compartment and insert 3x AA Batteries as illustrated overleaf
3. Connect monitor to Sensor Mat as illustrated (when the monitor stops beeping and the lights go out, the system is ready for testing.)
4. Position the Sensor in the most appropriate position (at the Carers discretion)
5. Attach the Ramblegard Monitor adjacent to or attached to the chair
6. TEST the system by pressing down on the pad. The LED should flash several times, and the monitor should beep twice to indicate "system in use". Release pressure from the pad and the monitor should alarm.
7. To reset the system either:
  - Press the reset button (or use optional carer key) to silence and set the system on standby.
  - or reapply pressure to the pad to silence and rearm the system.



full details overleaf

## Product Information

Product Code: SG1



**Product Name**  
Wired Seatgard

### System Contents

- 1x Ramblegard Corded Monitor - RG1
- 1x Ramblegard Corded Seatmat - RG10



Ramblegard have programmed the components included as a plug and play system.

If individual components need reprogrammed or replaced please refer to the relevant specification sheets available to download at [www.ramblegard.com](http://www.ramblegard.com)

**QUICKSTART INSTRUCTIONS** see overleaf

*Carers should also consider the following points when determining the most appropriate position for the the Ramblegard Sensor Pads & Monitors in care environments:*


#### Seat Sensor Pad Placement

- Position pad centred on the chair seat.

Monitor Placement - Attach using the quick release strap supplied or attach to wall or vertical surface using integral keyhole slots.

- The monitor can be placed within range of the connecting cable.
- The monitor should be placed where it cannot be tampered with by patients especially in dementia & challenging environment situations.

Back of Monitor



Disclaimer: Whilst every effort has been made to ensure the safety of this product, we cannot be held responsible for any damage or injury sustained when the equipment is in use. It is the responsibility of the user to install the system and ensure it is tested and operating correctly. The system is not designed to replace good care giving practice.