

Product Name: **Wireless Under-Mattress Bed Sensor Pad**

Product Code: **RG39**

Description: The Under-Mattress Bed Sensor Pad works with all Ramblegard wireless monitoring devices for both standalone and nurse call. This can be placed beneath the majority of dynamic and static mattresses. Triggers alarm when patient exits bed.

**Features:**

- Triggers alarm when pressure is released
- Eyelets for securing to a bed frame. Works under spring, static and dynamic mattresses 5-7 inches thick
- Durable and long-lasting
- Anti-bacterial PVC
- Wireless link with monitor or nurse call Jack
- Care instructions written on pad

**Works with:**

- RG2 Companion Monitor
- RG43 Wireless Display Panel
- RG41 Wireless Jill Pager
- RG5 Wireless Jack adaptor

**Technical Specifications**

<b>Connection</b>	4-Pin 2.5mm Transmitter Connector
<b>Material</b>	PVC & PET
<b>Dimensions</b>	76 x 25 cm
<b>Net Weight</b>	635g
<b>Operating Temp</b>	5°C to 40°C
<b>Operating Weight</b>	20 – 200 kg
<b>Internal Resistance</b>	<1200Ω
<b>Recommended Cleaning</b>	1:100 Solution of Household Bleach (5.25-61.5% Sodium Hypochlorite)

**Standards:**

EMC: EN55022:2010  
 Social Alarms: EN50130-4:1995  
 Manufacture: ISO9001:2003, ISO13485:2008

